

Anti-Cetn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41208

Item Number: CM0020086



Introduction

Rabbit polyclonal antibody against CETN2, validated for WB, IF/ICC, ELISA in human, mouse, rat. Involved in centrosome function, DNA repair, and mRNA export. Ideal for cell biology and genome stability studies.

[Learn More](#)

Parameter	Value
Product Name	CETN2 Rabbit pAb
Remarks/Alias	CALT; CEN2; CETN2
Species	Human
GeneID (human)	1069
GeneID	1069
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CETN2 (NP_004335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41208
Protein Weight	20kDa
Shipping	Ice bag