

Anti-Cetn3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15182

Item Number: CM0025730

Introduction

CETN3 rabbit polyclonal antibody validated for western blot, immunohistochemistry, immunofluorescence, and ELISA. Reacts with human, mouse, rat. Targets centrin-3, a centrosomal and nuclear pore component involved in mRNA export.

[Learn More](#)

Parameter	Value
Product Name	CETN3 Rabbit pAb
Remarks/Alias	CEN3; CDC31; CETN3
Species	Human
Gene ID (Human)	1070
Gene ID	1070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human CETN3 (NP_004356.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15182
Protein Weight	20kDa
Shipping	Ice bag