

Anti-Cetp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11597

Item Number: CM0014046



Introduction

Rabbit polyclonal antibody against human CETP for Western blot (WB), immunohistochemistry (IHC-P), and ELISA applications. Ideal for research on cholesterol metabolism, lipid transport, and atherosclerosis.

[Learn More](#)

Parameter	Value
Product Name	CETP Rabbit pAb
Remarks/Alias	BPIFF; HDLCQ10; CETP
Species	Human
GeneID (Human)	1071
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-493 of human CETP (NP_000069.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11597
Protein Weight	55kDa
Shipping	Ice bag