

Anti-Cfap20 (C16Orf80) Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6A4

Item Number: CM0026789



Introduction

Rabbit polyclonal antibody targeting human C16orf80/CFAP20, validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studying ciliary motility and axonemal organization.

[Learn More](#)

Parameter	Value
Product Name	C16orf80 Rabbit pAb
Remarks/Alias	GTL3; BUG22; EVORF; fSAP23; C16orf80
Species	Human
GeneID (human)	29105
GeneID	29105
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 94-193 of human C16orf80 (NP_037374.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6A4
Protein Weight	23kDa
Shipping	Ice bag