


Anti-Cfap44 Polyclonal Antibody For Wb, Elisa - Q96Mt7

Item Number: CM0013828



Introduction

Rabbit polyclonal antibody against human CFAP44 (WDR52) for Western blot and ELISA applications. Cross-reacts with mouse. Suitable for sperm flagellum axoneme organization studies.

[Learn More](#)

Parameter	Value
Product Name	CFAP44 Rabbit pAb
Remarks/Alias	WDR52; SPGF20; CFAP44
Species	Human
GeneID (Human)	55779
GeneID	55779
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-370 of human CFAP44 (NP_001157968.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96MT7
Protein Weight	214kDa
Shipping	Ice bag