

Anti-Cfap53 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96M91

Item Number: CM0021975



Introduction

High-quality rabbit polyclonal antibody targeting CFAP53 (CCDC11/HTX6), suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat CFAP53. Essential for studying motile cilia and flagella biogenesis.

[Learn More](#)

Parameter	Value
Product Name	CFAP53 Rabbit pAb
Remarks/Alias	HTX6; CCDC11; CFAP53
Species	Human
Gene ID (Human)	220136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CFAP53 (NP_659457.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96M91
Protein Weight	62kDa
Shipping	Ice bag