


## Anti-Cfd Rabbit Monoclonal Antibody For Wb, Elisa - P00746

Item Number: CM0003931



### Introduction

Rabbit monoclonal antibody against human complement factor D (Adipsin). Validated for WB and ELISA applications. Ideal for innate immunity, complement pathway studies, and diagnostic raw material sourcing.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | CFD Rabbit mAb                                                                                             |
| Remarks/Alias    | DF; ADN; PFD; ADIPSIN; CFD                                                                                 |
| Species          | Human                                                                                                      |
| GeneID (Human)   | 1675                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-253 of human CFD (P00746) |
| Source           | Rabbit                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                      |
| Application      | WB, ELISA                                                                                                  |
| Cross Reactivity | Human                                                                                                      |
| SWISS            | P00746                                                                                                     |
| Protein Weight   | 27kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                    |