


Anti-Cfd Rabbit Polyclonal Antibody For Wb, Elisa - P00746

Item Number: CM0027608



Introduction

Rabbit polyclonal antibody against human CFD (Adipsin, Complement Factor D). Validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for studying the alternative complement pathway and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	CFD Rabbit pAb
Remarks/Alias	DF; ADN; PFD; ADIPSIN; CFD
Species	Human
GeneID (human)	1675
GeneID	1675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-253 of human CFD (NP_001919.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00746
Protein Weight	27kDa
Shipping	Ice bag