

Anti-Cfdp1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uee9

Item Number: CM0026049

Introduction

Anti-CFDP1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets heterochromatin-stabilizing protein involved in centromeric stability and craniofacial development. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	CFDP1 Rabbit pAb
Remarks/Alias	p97; BCNT; CP27; SWC5; Yeti; CENP-29; BUCENTAUR; CFDP1
Species	Human
GeneID	10428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human CFDP1 (NP_006315.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UEE9
Protein Weight	34kDa
Shipping	Ice bag