


Anti-Cfhr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q03591

Item Number: CM0024453



Introduction

Rabbit polyclonal antibody targeting human CFHR1, suitable for Western blot and ELISA with cross-reactivity to mouse and rat. Recognizes a ~38kDa secreted glycoprotein involved in complement regulation and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	CFHR1 Rabbit pAb
Remarks/Alias	H36; CFHL; FHR1; HFL1; HFL2; CFHL1; FHL-1; FHR-1; H36-1; H36-2; CFHL1P; CFHR1P; CFHR1
Species	Human
Gene ID (Human)	3078
Gene ID	3078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 151-330 of human CFHR1 (NP_002104.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03591
Protein Weight	38kDa
Shipping	Ice bag