

Anti-Cfhr3 Polyclonal Antibody For Wb, If-P, Elisa - Q02985

Item Number: CM0026741

Introduction

This CFHR3 Rabbit Polyclonal Antibody is validated for Western blot, IF-P, and ELISA in human, mouse, and rat. Ideal for studying complement regulation and hemolytic uremic syndrome.

[Learn More](#)

Parameter	Specification
Product Name	CFHR3 Rabbit pAb
Remarks/Alias	FHR3; HLF4; CFHL3; FHR-3; DOWN16; CFHR3
Species	Human
GeneID (Human)	10878
GeneID	10878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-200 of human CFHR3 (NP_066303.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02985
Protein Weight	37kDa
Shipping	Ice bag