

Anti-Cfhr5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bxr6

Item Number: CM0017322

Introduction

High-quality rabbit polyclonal antibody targeting human CFHR5 (FHR5), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for complement regulation and hemolytic uremic syndrome research.

[Learn More](#)

Parameter	Value
Product Name	CFHR5 Rabbit pAb
Remarks/Alias	FHR5; CFHL5; FHR-5; CFHR5D; CFHR5
Species	Human
Gene ID (Human)	81494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-569 of human CFHR5 (NP_110414.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXR6
Protein Weight	64kDa
Shipping	Ice bag