


Anti-Cfi Polyclonal Antibody For Wb And Elisa - P05156

Item Number: CM0026615



Introduction

Anti-CFI rabbit polyclonal antibody suitable for WB and ELISA applications, reactive with human and mouse CFI. Ideal for complement pathway studies, immune regulation research, and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CFI Rabbit pAb
Remarks/Alias	FI; IF; KAF; AHUS3; ARMD13; C3BINA; C3b-INA; CFI
Species	Human
GeneID (Human)	3426
GeneID	3426
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-300 of human CFI (NP_000195.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05156
Protein Weight	66kDa
Shipping	Ice bag