


Anti-Cfl1 Rabbit Polyclonal Antibody For Wb And Elisa - P23528

Item Number: CM0020897



Introduction

CFL1 Rabbit pAb is a rabbit polyclonal antibody against human cofilin-1, validated for WB and ELISA. It detects human, mouse, and rat CFL1 (predicted 19kDa), ideal for studying actin dynamics, cell migration, and cytoskeleton organization.

[Learn More](#)

Parameter	Value
Product Name	CFL1 Rabbit pAb
Remarks/Alias	CFL; cofilin; HEL-S-15; CFL1
Species	Human
GeneID (human)	1072
GeneID	1072
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CFL1 (NP_005498.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23528
Protein Weight	19kDa
Shipping	Ice bag