

Anti-Cfl2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y281

Item Number: CM0026036



Introduction

Cofilin-2 rabbit polyclonal antibody validated for Western blotting, immunofluorescence, and ELISA. Reacts with human, mouse, and rat CFL2. Ideal for muscle actin dynamics and nemaline myopathy research.

[Learn More](#)

Parameter	Value
Product Name	CFL2 Rabbit pAb
Remarks/Alias	NEM7; CFL2
Species	Human
GeneID (Human)	1073
GeneID	1073
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human CFL2 (NP_068733.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y281
Protein Weight	19kDa
Shipping	Ice bag