


Anti-Cflip Polyclonal Antibody For Wb, Elisa - 015519

Item Number: CM0023326



Introduction

Rabbit polyclonal antibody against human cFLIP (CASP8 and FADD-like apoptosis regulator) for Western blot, ELISA. Cross-reacts with human, mouse, rat. Suitable for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	cFLIP Rabbit pAb
Remarks/Alias	CASH; FLIP; MRIT; CLARP; FLAME; cFLIP; Casper; FLAME1; c-FLIP; FLAME-1; I-FLICE; c-FLIPL; c-FLIPR; c-FLIPS; CASP8AP1
Species	Human
GeneID	8837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-240 of human cFLIP (NP_003870.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	015519
Protein Weight	55kDa
Shipping	Ice bag