


Anti-Cfp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27918

Item Number: CM0027311



Introduction

Rabbit polyclonal antibody targeting human Properdin (CFP), validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat. Recombinant immunogen. For complement alternate pathway research.

[Learn More](#)

Parameter	Value
Product Name	CFP Rabbit pAb
Remarks/Alias	BFD; PFC; PFD; PROPERDIN; CFP
Species	Human
Gene ID (Human)	5199
Gene ID	5199
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 190-469 of human CFP (NP_001138724.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27918
Protein Weight	51kDa
Shipping	Ice bag