


Anti-Cftr Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13569

Item Number: CM0026556



Introduction

Rabbit polyclonal antibody targeting human CFTR, validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, rat CFTR (~168 kDa). Ideal for ion channel and cystic fibrosis research.

[Learn More](#)

Parameter	Value
Product Name	CFTR Rabbit pAb
Remarks/Alias	CF; MRP7; ABC35; ABCC7; CFTR/MRP; TNR-CFTR; dJ760C5.1; CFTR
Species	Human
Gene ID	1080
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-840 of human CFTR (NP_000483.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13569
Protein Weight	168kDa
Shipping	Ice bag