

Anti-Cgas Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N884

Item Number: CM0030589



Introduction

Rabbit polyclonal antibody against human cGAS (cyclic GMP-AMP synthase; MB21D1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for innate immunity and DNA sensing research.

[Learn More](#)

Parameter	Value
Product Name	cGAS Rabbit pAb
Remarks/Alias	MB21D1; h-cGAS; C6orf150; cGAS
Species	Human
GeneID (Human)	115004
GeneID	115004
Immunogen	synthetic peptide Recombinant fusion protein containing a sequence corresponding to amino acids 350-522 of human cGAS (NP_612450.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N884
Protein Weight	59 kDa
Shipping	Ice bag