

## Anti-Cgas Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N884

Item Number: CM0008929



### Introduction

Rabbit monoclonal antibody targeting human cGAS (MB21D1, C6orf150), suitable for WB, IHC-P, IF/ICC, ELISA. Recognizes a band of ~59 kDa. Key reagent for innate immunity and DNA sensing studies.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | cGAS Rabbit mAb                                                                                        |
| Remarks/Alias    | MB21D1; h-cGAS; C6orf150                                                                               |
| Species          | Human                                                                                                  |
| Gene ID (Human)  | 115004                                                                                                 |
| Gene ID          | 115004                                                                                                 |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 430-522 of human cGAS (NP_612450.2) |
| Source           | Rabbit                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                  |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                               |
| Cross Reactivity | Human                                                                                                  |
| SWISS            | Q8N884                                                                                                 |
| Protein Weight   | 59 kDa                                                                                                 |
| Shipping         | Ice bag                                                                                                |