

Anti-Cgb Rabbit Polyclonal Antibody For Wb, Elisa - P0Dn86

Item Number: CM0017215



Introduction

High-quality rabbit polyclonal antibody against human CGB (CGB3/CGB5/CGB8), validated for Western blot and ELISA. Detects Human, Mouse, Rat. Ideal for pregnancy research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CGB Rabbit pAb
Remarks/Alias	CGB; LHB; CGB5; CGB7; CGB8; hCGB
Species	Human
GeneID (Human)	1082
GeneID	1082
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DN86
Protein Weight	18kDa
Shipping	Ice bag