

## Anti-Cgb5 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dn86

Item Number: CM0029195



### Introduction

Rabbit polyclonal antibody against human CGB5, suitable for WB and ELISA. Validated for human samples, targeting the beta subunit of hCG essential in pregnancy and maternal adaptation.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CGB5 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | CGB; HCG; hCGB; CGB5                                                                                                                  |
| Species          | Human                                                                                                                                 |
| GeneID           | 93659                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-165 of human CGB5 (NP_000728.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P0DN86                                                                                                                                |
| Protein Weight   | 18kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |