

Anti-Cgb7 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dn87

Item Number: CM0029419

Introduction

CGB7 Rabbit pAb is a polyclonal antibody raised in rabbit against human CGB7. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Suitable for pregnancy research and hCG studies.

[Learn More](#)

Parameter	Value
Product Name	CGB7 Rabbit pAb
Remarks/Alias	CGB6; CG-beta-a; CGB7
Species	Human
Gene ID (Human)	94027
Gene ID	94027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-165 of human CGB7 (NP_149133.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DN87
Protein Weight	18kDa
Shipping	Ice bag