

Anti-Cgb7 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dn87

Item Number: CM0015221

Introduction

High-quality rabbit polyclonal antibody against human Choriogonadotropin subunit beta 7 (CGB7). Validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, rat. Ideal for reproductive biology and pregnancy research.

[Learn More](#)

Parameter	Value
Product Name	CGB7 Rabbit pAb
Remarks/Alias	CGB6; CG-beta-a; CGB7
Species	Human
GeneID	94027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-165 of human CGB7 (NP_149133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DN87
Protein Weight	18kDa
Shipping	Ice bag