

Anti-Cggbp1 Polyclonal Antibody For Wb, Elisa - Q9Ufw8

Item Number: CM0029616

Introduction

CGGBP1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human CGGBP1, a DNA-binding protein that associates with CGG repeats in the FMR1 promoter. 19 kDa.

[Learn More](#)

Parameter	Value
Product Name	CGGBP1 Rabbit pAb
Remarks/Alias	CGGBP; p20-CGGBP; CGGBP1
Species	Human
Gene ID (Human)	8545
Gene ID	8545
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human CGGBP1 (NP_003654.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UFW8
Protein Weight	19kDa
Shipping	Ice bag