

Anti-Cgref1 Rabbit Polyclonal Antibody For Wb, Elisa - Q99674

Item Number: CM0015680

Introduction

Rabbit polyclonal antibody targeting CGREF1, a cell growth regulator and adhesion protein. Suitable for WB and ELISA, cross-reacts with mouse and rat. Useful for cancer and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	CGREF1 Rabbit pAb
Remarks/Alias	CGR11; CGREF1
Species	Human
GeneID (human)	10669
GeneID	10669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human CGREF1 (NP_006560.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99674
Protein Weight	33kDa
Shipping	Ice bag