

Anti-Cgrp Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P01258 / P06881

Item Number: CM0030279



Introduction

Rabbit polyclonal antibody against human CGRP (calcitonin gene-related peptide), suitable for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for calcium metabolism and neuropeptide research.

[Learn More](#)

Parameter	Value
Product Name	CGRP Rabbit pAb
Remarks/Alias	CT, KC, PCT, CGRP, CALC1, CGRP1, CGRP-I, CGRP-alpha
Species	Human
GeneID (Human)	796
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 49-128 of human CGRP (NP_001029125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01258/P06881
Protein Weight	13kDa/15kDa
Shipping	Ice bag