

Anti-Chac1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Bux1

Item Number: CM0016797

Introduction

Rabbit polyclonal antibody targeting human CHAC1. Suitable for WB, IF-P, IHC-P, and ELISA. Detects CHAC1, a pro-apoptotic enzyme implicated in unfolded protein response and Notch pathway regulation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CHAC1 Rabbit pAb
Remarks/Alias	CHAC1
Species	Human
Gene ID (Human)	79094
Gene ID	79094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human CHAC1 (NP_077016.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUX1
Protein Weight	24kDa
Shipping	Ice bag