

Anti-Chad Rabbit Polyclonal Antibody For Wb, Elisa - O15335

Item Number: CM0018674



Introduction

Rabbit polyclonal antibody against human chondroadherin (CHAD) protein, suitable for western blot (WB) and ELISA. Useful for investigating cartilage ECM biology and chondrocyte adhesion mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CHAD Rabbit pAb
Remarks/Alias	SLRR4A; CHAD
Species	Human
GeneID (Human)	1101
GeneID	1101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human CHAD (NP_001258.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15335
Protein Weight	40kDa
Shipping	Ice bag