

Anti-Chat Polyclonal Antibody For Wb, Ihc, If, Elisa - P28329

Item Number: CM0018231



Introduction

High-quality rabbit polyclonal antibody against CHAT (choline O-acetyltransferase), validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat CHAT with specificity. Ideal for neuroscience and neuromuscular research.

[Learn More](#)

Parameter	Value
Product Name	CHAT Rabbit pAb
Remarks/Alias	CMS6; CMS1A; CMS1A2; CHOACTASE; CHAT
Species	Human
Gene ID (Human)	1103
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 331-630 of human CHAT (NP_065574.4))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28329
Protein Weight	83kDa
Shipping	Ice bag