

Anti-Chat Rabbit Monoclonal Antibody For Wb, Elisa - P28329

Item Number: CM0002727



Introduction

Recombinant rabbit monoclonal antibody targeting human choline O-acetyltransferase (ChAT). Validated for WB and ELISA with human, mouse, and rat reactivity. Ideal for cholinergic research.

[Learn More](#)

Parameter	Value
Product Name	ChAT Rabbit mAb
Remarks/Alias	CMS6; CMS1A; CMS1A2; CHOACTASE; CHAT
Species	Human
GeneID (Human)	1103
GeneID	1103
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28329
Protein Weight	83kDa
Shipping	Ice bag