


Anti-Chat Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P28329

Item Number: CM0022648



Introduction

Rabbit polyclonal antibody targeting human CHAT (Choline O-acetyltransferase), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. ~83kDa. Ideal for cholinergic synapse and neurotransmitter studies.

[Learn More](#)

Parameter	Value
Product Name	CHAT Rabbit pAb
Remarks/Alias	CMS6; CMS1A; CMS1A2; CHOACTASE; CHAT
Species	Human
GeneID (Human)	1103
GeneID	1103
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CHAT (NP_065574.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28329
Protein Weight	83kDa
Shipping	Ice bag