

Anti-Chchd3 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Nx63

Item Number: CM0028786

Introduction

Rabbit polyclonal antibody against human CHCHD3 (MIC19/MINOS3). Validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial research and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CHCHD3 Rabbit pAb
Remarks/Alias	Mic19; MINOS3; MICOS19; PPP1R22; CHCHD3
Species	Human
Gene ID (Human)	54927
Gene ID	54927
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX63
Protein Weight	26kDa
Shipping	Ice bag