

Anti-Chchd3 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q9Nx63

Item Number: CM0022857

Introduction

KO-validated rabbit polyclonal antibody against human CHCHD3 (Mic19), suitable for WB and ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial MICOS complex and cristae morphology studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CHCHD3 Rabbit pAb
Remarks/Alias	Mic19; MINOS3; MICOS19; PPP1R22; D3
Species	Human
GeneID (Human)	54927
GeneID	54927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human CHCHD3 (NP_060282.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX63
Protein Weight	26kDa
Shipping	Ice bag