

Anti-Chchd4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N4Q1

Item Number: CM0018901



Introduction

Rabbit polyclonal antibody against human CHCHD4/MIA40, validated for WB, IHC-P, and ELISA. Ideal for studying mitochondrial intermembrane space import and redox-sensitive protein folding; cross-reacts with rat.

[Learn More](#)

Parameter	Value
Product Name	CHCHD4 Rabbit pAb
Remarks/Alias	MIA40; TIMM40; CHCHD4
Species	Human
GeneID (Human)	131474
GeneID	131474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human CHCHD4 (NP_001091972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8N4Q1
Protein Weight	16kDa
Shipping	Ice bag