


Anti-Chchd4 Polyclonal Antibody (Rabbit Pab) For Wb, Ihc-P, Elisa - Q8N4Q1

Item Number: CM0014123



Introduction

High-quality rabbit polyclonal antibody against human CHCHD4 (MIA40), validated for WB, IHC-P, and ELISA. Suitable for studying mitochondrial protein import and respiratory chain assembly. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	CHCHD4 Rabbit pAb
Remarks/Alias	MIA40; TIMM40; CHCHD4
Species	Human
GeneID (Human)	131474
GeneID	131474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human CHCHD4 (NP_001091972.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4Q1
Protein Weight	16kDa
Shipping	Ice bag