


Anti-Chchd4 Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q8N4Q1

Item Number: CM0012237



Introduction

Rabbit monoclonal antibody against CHCHD4 (MIA40), validated for WB, IHC, IF, and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CHCHD4 Rabbit mAb
Remarks/Alias	MIA40; TIMM40
Species	Human
Gene ID (Human)	131474
Gene ID	131474
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human CHCHD4 (NP_001091972.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4Q1
Protein Weight	16kDa
Shipping	Ice bag