

Anti-Chchd6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Brq6

Item Number: CM0017403

Introduction

Rabbit polyclonal antibody against human CHCHD6 (MICOS complex subunit MIC25), suitable for WB and ELISA. Targets a mitochondrial inner membrane protein crucial for crista junction architecture. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CHCHD6 Rabbit pAb
Remarks/Alias	CHCM1; Mic25; MICOS25; PPP1R23; CHCHD6
Species	Human
Gene ID (Human)	84303
Gene ID	84303
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-235 of human CHCHD6 (NP_115719.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BRQ6
Protein Weight	26kDa
Shipping	Ice bag