

Anti-Chchd7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Buk0

Item Number: CM0023515



Introduction

High-quality rabbit polyclonal antibody targeting human CHCHD7 (COX23). Validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat. Ideal for mitochondrial research and chromosomal rearrangement studies.

[Learn More](#)

Parameter	Value
Product Name	CHCHD7 Rabbit pAb
Remarks/Alias	COX23; CHCHD7
Species	Human
Gene ID	79145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human CHCHD7 (NP_077276.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUK0
Protein Weight	10kDa
Shipping	Ice bag