

Anti-Chd1 Rabbit Monoclonal Antibody For Wb, Elisa - O14646

Item Number: CM0005326



Introduction

Anti-CHD1 Rabbit Monoclonal Antibody validated for WB and ELISA. Targets human CHD1, an ATP-dependent chromatin remodeler involved in transcription and splicing. Ideal for epigenetic and developmental research.

[Learn More](#)

Parameter	Value
Product Name	CHD1 Rabbit mAb
Remarks/Alias	CHD-1; PILBOS; CHD1
Species	human
GeneID	1105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1510-1710 of human CHD1 (NP_001261.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14646
Protein Weight	197 kDa
Shipping	Ice bag