

Anti-Chd2 Polyclonal Antibody For Wb, Elisa - O14647

Item Number: CM0024747

Introduction

Validate CHD2 chromatin remodeling function with this rabbit polyclonal antibody, suitable for WB and ELISA in human, mouse, and rat samples. Ideal for epigenetics and myogenesis research.

[Learn More](#)

Parameter	Value
Product Name	CHD2 Rabbit pAb
Remarks/Alias	EEOC; DEE94; CHD2
Species	Human
Gene ID (human)	1106
Gene ID	1106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human CHD2 (NP_001036037.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14647
Protein Weight	211kDa
Shipping	Ice bag