

Anti-Chd2 Polyclonal Antibody For Wb, Elisa - O14647

Item Number: CM0014021

Introduction

CHD2 Rabbit Polyclonal Antibody validated for WB and ELISA against human CHD2. Ideal for chromatin remodeling, myogenesis, and transcription regulation studies in human cell lines.

[Learn More](#)

Parameter	Value
Product Name	CHD2 Rabbit pAb
Remarks/Alias	EEOC, DEE94, CHD2
Species	Human
Gene ID (Human)	1106
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1729-1828 of human CHD2 (NP_001262.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14647
Protein Weight	211kDa
Shipping	Ice bag