

Anti-Chd3 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q12873

Item Number: CM0007931



Introduction

Anti-CHD3 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Detects human, mouse, rat CHD3, a 227kDa NuRD complex chromatin remodeler involved in transcriptional repression. Ideal for epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	CHD3 Rabbit mAb
Remarks/Alias	ZFH; Mi-2a; SNIBCP5; Mi2-ALPHA; CHD3
Species	Human
GenelD (Human)	1107
GenelD	1107
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12873
Protein Weight	227kDa
Shipping	Ice bag