


Anti-Chd4 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14839

Item Number: CM0007135



Introduction

Rabbit monoclonal antibody against human CHD4, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Target is an ATP-dependent chromatin remodeler in the NuRD complex, involved in transcription regulation and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	CHD4 Rabbit mAb
Remarks/Alias	CHD-4; Mi-2b; SIHIWES; Mi2-BETA; CHD4
Species	Human
Gene ID	1108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1520-1700 of human CHD4. (NP_001264.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14839
Protein Weight	218kDa
Shipping	Ice bag