

Anti-Chd9 Polyclonal Antibody For Wb And Elisa - Q3L8U1

Item Number: CM0018320

Introduction

CHD9 Rabbit pAb - validated for WB and ELISA. Detects human and mouse CHD9 (predicted 326 kDa), a chromatin remodeler involved in transcriptional regulation. Recombinant immunogen (aa 85-190). For research use.

[Learn More](#)

Parameter	Value
Product Name	CHD9 Rabbit pAb
Remarks/Alias	AD013; CHD-9; CReMM; KISH2; PRIC320; CHD9
Species	Human
Gene ID (Human)	80205
Gene ID	80205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 85-190 of human CHD9 (NP_079410.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3L8U1
Protein Weight	326kDa
Shipping	Ice bag