

Anti-Chdh Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ne62

Item Number: CM0019288

Introduction

Rabbit polyclonal antibody targeting human CHDH (choline dehydrogenase, mitochondrial). Validated for WB, IHC-P, ELISA, and detects human, mouse, rat species. Suitable for studying choline metabolism, mitochondrial autophagy, and related pathways.

[Learn More](#)

Parameter	Value
Product Name	CHDH Rabbit pAb
Remarks/Alias	CHDH
Species	Human
GeneID (Human)	55349
GeneID	55349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-594 of human CHDH (NP_060867.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NE62
Protein Weight	65kDa
Shipping	Ice bag