

Anti-Chek2 Polyclonal Antibody For Wb, Elisa - Q9Z265

Item Number: CM0020193

Introduction

High-quality rabbit polyclonal antibody against human Chek2, validated for Western blotting and ELISA with human and mouse cross-reactivity. Ideal for DNA damage response and cell cycle checkpoint studies in cancer research.

[Learn More](#)

Parameter	Value
Product Name	Chek2 Rabbit pAb
Remarks/Alias	CHK2; Cds1; Rad53; HUCDS1; Chek2
Species	Human
GeneID (Human)	50883
GeneID	50883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human Chek2. (NP_057890.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Z265
Protein Weight	61kDa
Shipping	Ice bag