

Anti-Chemerin Rabbit Polyclonal Antibody For Wb, Elisa - Q9Dd06

Item Number: CM0030257

Introduction

Rabbit polyclonal antibody against mouse Chemerin (Rarres2), validated for Western blot and ELISA. Reacts with mouse and rat. Suitable for adipokine, inflammation, and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Chemerin Rabbit pAb
Remarks/Alias	0610007L05Rik
Species	Mouse
Gene ID (Human)	5919
Gene ID	71660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-162 of mouse Chemerin (NP_082128.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9DD06
Protein Weight	18kDa
Shipping	Ice bag