

Anti-Cherp Polyclonal Antibody For Wb, Elisa - Q8Iwx8

Item Number: CM0030596

Introduction

CHERP Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat CHERP protein. Ideal for studies on calcium homeostasis and endoplasmic reticulum function.

[Learn More](#)

Parameter	Value
Product Name	CHERP Rabbit pAb
Remarks/Alias	SRA1; DAN16; SCAF6; CHERP
Species	Human
Gene ID (Human)	10523
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 650-750 of human CHERP (NP_006378.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWX8
Protein Weight	104kDa
Shipping	Ice bag