

## Anti-Chga Polyclonal Antibody For Wb, Elisa - P10645

Item Number: CM0020780



### Introduction

Rabbit polyclonal CHGA antibody validated for WB and ELISA. Reacts with human and mouse CHGA (51 kDa), a granin family protein regulating hormone secretion and innate immunity. Suitable for neuroendocrine and cardiovascular research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CHGA Rabbit pAb                                                                                                    |
| Remarks/Alias    | CGA; PHES; PHES; CHGA                                                                                              |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 1113                                                                                                               |
| Gene ID          | 1113                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 230-380 of human CHGA (NP_001266.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | P10645                                                                                                             |
| Protein Weight   | 51kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |